

Standard Analysis Certificate

PRODUCT :p-Aminobenzoic acid ERBApharm-According to pharmacopoeia: USP
CODE :391800
METHOD :6098

TEST	U.M.	SPECIFICATION
Description	-	Yellowish powder
Identification	-	Positive
Ordinary impurities	USP-NF	Conform
Melting point	°C	186 ÷ 189
Loss on drying	%	<=0.2
Sulphated ash	%	<=0.1
Volat. diazotizable sub	ppm	<=20
Heavy metals (Pb)	ppm	<=20
Assay (acidimetric)	% on dry	98.5 ÷ 100.5

Our products should be used in compliance with the current legislation, raw material for pharmaceutical uses included.

Date :14/07/2018

ЛАБОРАТОРКА
с нами делают открытия

QUALITY CONTROL RESPONSIBLE
B. COULANGE (VDR)

ООО "Лабораторка" +7 916 869-98-67, contact@laboratorka.com, web: www.laboratorka.com