

Standard Analysis Certificate

PRODUCT :Aluminum potassium sulfate dodecahydrate ERBApharm-According to pharmacopoeia: Ph.Eur.-USP-FU-Ph.Franc.-BP
CODE :312400
METHOD :7357

TEST	U.M.	SPECIFICATION
Description	-	White crystalline powder
Identification	-	Positive
Appearance of solution	Ph.Eur.	Conform
pH 10% at 25° C	-	3.0 ÷ 3.5
Loss on drying	%	43.0 ÷ 46.0
Ammonium	%	<= 0.2
Heavy metals (Pb)	ppm	<= 20
Fe	ppm	<= 100
Assay (complexometric)	%	99.0 ÷ 100.5

Our products should be used in compliance with the current legislation, raw material for pharmaceutical uses included.

Date :14/07/2018