

Standard Analysis Certificate

PRODUCT :ETHYLENEDIAMINETETRAACETIC ACID DISODIUM SALT ERBApharm-According to pharmacopoeia: Ph.Eur.-USP

CODE :303220

METHOD :7425

TEST	U.M.	SPECIFICATION
Description	-	White powder
Identification	-	Positive
Appearance of solution	Ph.Eur.	Conform
Ca	USP-NF	Conform
pH sol. 5% at 25° C	-	4.0 ÷ 5.5
Loss on drying	%	8.7 ÷ 11.4
Ac. nitrilotriacetic (Impurezza A)	%	<= 0.1
Heavy metals (Pb)	ppm	<= 20
Fe	ppm	<= 80
Assay (complexometric)	% on dry	99.0 ÷ 101.0
Assay (complexometric)	% t.q.	98.5 ÷ 101.0
Origin (BSE/TSE)	-	Synthesis
Residual solvents (Current ICH)	-	Conform

Our products should be used in compliance with the current legislation, raw material for pharmaceutical uses included.

Date

:06/01/2018

QUALITY CONTROL RESPONSIBLE
B. COULANGE (VDR)